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Transformation of *E. coli* with plasmid DNA → Extraction of the plasmid

Day 1: *E. coli* Transformation with Plasmid DNA (evening onwards)

Materials:

- Chemical competent cells
- Plasmid DNA or ligation product
- Agar plates with an antibiotic
- SOC medium or LB medium
- ColiRoller™ plating beads or Bacterial spreaders

- 1) Preheat the heat block to 42 °C.
- 2) Retrieve tubes containing competent cells (10 µL to 100 µL) from –80 °C (Deep freezer) and allow them to thaw on ice. *If the volume is around 100 µL, you can briefly thaw them at room temperature before placing them on ice until partially melted.
- 3) Once the competent cells have thawed, add the DNA and incubate on ice for 30 minutes. *For larger volumes of competent cell solution (≥ 50 µL), gently pipette to mix the cells and DNA several times. Conversely, if the volume is 10 µL, mixing is not necessary. For cyclic plasmid DNA transformation, a 15-minute incubation on ice may suffice.
- 4) Perform a heat shock in the preheated heat block at 42 °C for 40–60 seconds.
- 5) Immediately chill the cell on ice and let it sit for approximately 2 minutes.
- 6) Add an appropriate volume of SOC or LB medium, gently pipette, and spread all or part of the cell suspension onto the agar plate. *For Ampicillin resistance, immediately spread *E. coli* onto an ampicillin agar plate. For Kanamycin resistance, preincubate in antibiotic-free medium at 37 °C for 1 hour before spreading onto the plate (critical step!). This arises from differences in antibiotic action mechanisms (please investigate the reason by yourself).
- 7) Incubate the plates with the agar side up at 37 °C overnight. *After 12–16 hours, if no colonies are visible, consider the transformation unsuccessful. If there are even a few and small colonies, continue incubation for a bit longer.

Note1:

Transformation efficiency varies based on whether the competent cells are homemade or purchased and whether the DNA is cyclic plasmid DNA or a ligation product. Typically, 0.1–20 pg of DNA is used. For plasmid DNA transformation, 1–10 pg is sufficient, while for ligation products, 20–50 ng of DNA (amount of linearized plasmid DNA) is recommended.

Note2:

Handle *E. coli* very gently and avoid warming them until the heat shock step.

Note3:

Different strains of *E. coli* serve different purposes, such as amplifying and extracting plasmids or expressing and purifying proteins. Choose the appropriate strain accordingly. Common strains include "JM109," "DH5 α ," and "XL1-Blue" for plasmid preparation and BL21 series for protein expression.

Day 2: Colony picking (morning & evening)**Materials:**

- LB medium with antibiotics
 - Spitz tubes
- 1) Retrieve the agar plates inoculated with *E. coli* from the incubator (after 12 to 16 hours) and place them at room temperature or in a refrigerator. *During hot summer days, when the colonies grow densely, it's advisable to store them in the refrigerator to prevent clustering and sticking together. If colony picking is not immediate, refrigeration is recommended too. Cover the plates with plastic wrap to prevent them from drying out.
 - 2) Dispense 2–3 ml of LB medium with antibiotics into a Spitz tube.
 - 3) Using a sterile toothpick or a yellow pipette tip (200 μ L), pick the colonies from 1) and inoculate them into the medium in 2). Touch a single, well-isolated colony on the plate. *For very small colonies, a toothpick may be preferable as its rough surface reduces the risk of dislodging the colony. However, pipette tips can be also used. When picking a colony, you don't have to get it all off.
 - 4) Incubate overnight at 37°C with shaking.

Note1:

Pick up multiple colonies per plasmid DNA (or ligation product). When attempting to create a DNA construct containing the gene of interest (insert), it's common to encounter clones lacking the insert or containing mutations. By selecting multiple clones, you can increase the chances of obtaining the desired construct.

Day 3: Plasmid DNA Extraction (morning)**Materials:**

- Plasmid extraction kit (available from QIAGEN, Promega, MACHEREY-NAGEL, and other companies)
- 1) Centrifuge the *E. coli* suspension at 4 °C for about 10 minutes, then discard the supernatant by decantation.
 - 2) Follow the protocol provided with the kit. *The procedure is outlined as follows:
E. coli Suspension in Buffer (with RNase A) \rightarrow *E. coli* lysis in alkaline solution (NaOH-SDS), DNA

denaturation (to single strand) → Neutralization (KOAc) *where genomic DNA precipitates along with protein and SDS, while the plasmid DNA returns to be double-stranded and stay in the soluble fraction → Column washing → Elution of the plasmid DNA.

Note1:

Do not vortex after adding the alkaline solution, as it will shear the genomic DNA!